

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
Due Date:	

Scan Page 1 & 2 and email to: hs-cppn@ucdavis.edu Preceptor or other verified personnel date and initial completed skill.

Skill/Learning Not all skills are applicable to all Nursing areas – if not applicable mark as N/A	Completed Online Module	Date Completed (or N/A)	Verifier Initials
AdenoPlus by RPS Checklist 2024	DAHS-NSCAPRPS24-POCT		
Bacterial Vaginosis (OSOM BV Blue) Test Checklist 2024	DAHS-NSCBVRT24-POCT		
CoaguChek Checklist 2024	DAHS-NSCCOAGU24-POCT		
Cobas® Liat® System Checklist 2024	DAHS-NSCCOBAS24-POCT		
ColoScreen Checklist 2024	DAHS-NSCCOLO24-POCT		
Gastrocult Checklist 2024	DAHS-NSCGAST24-POCT		
HemoCue Hb 201 DM Checklist 2024	DAHS-NSCHEMOC24-POCT		
Hemoglobin A1c (HbA1c) Afinion Checklist 2024	DAHS-NSCHEMOG24-POCT		
Lactate Plus Checklist 2024	DAHS-NSCLPMSM24-POCT		
Medtronic HMS Plus	DAHS-NSCMEDHMSP24-POCT		
Nova StatStrip Checklist 2024	DAHS-NSCNOVA24-POCT		
pH Paper Checklist 2024	DAHS-NSCPHP24-POCT		
Urine Dipstick Clinitek Status+ Connect Checklist 2024	DAHS-NSCURID24-POCT		
Urine Drug Screen Alere iCup Checklist 2024	DAHS-NSCALERE24-POCT		
Urine Pregnancy Test Checklist 2024	DAHS-NSCURIP24-POCT		

Scan Page 1 and 2 and email to: hs-cppn@ucdavis.edu

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
Due Date:	

Scan Page 1 & 2 and email to: hs-cppn@ucdavis.edu Preceptor or other verified personnel date and initial completed skill.

SIGNATURE PAGE:

Signature and Printed Name of Verifier (preceptor or other verified personnel) who have initialed on this form:

Initial:	Print Name:	Signature:

PRECEPTEE STATEMENT AND SIGNATURE:

I have read and understand the appropriate UCDH Policies and Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills as noted, and I have the knowledge of the resources available to answer questions.

Printed Name	Signature

Scan Page 1 and 2 and email to: hs-cppn@ucdavis.edu

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:

These skills will be considered complete when all below performance criteria are completed.

AdenoPlus	# DAHS-NSCAPRPS24-POCT		
REFERENCES:		Date	Verifier Initials
1. UC Davis Health Policy POCT.57: AdenoPlus by RPS 2. RPS Package Insert AdenoPlus External Controls. Rapid Pathogen Screening Inc., Sarasota ,FL 34240			
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes sample collection and handling. <input type="checkbox"/> Demonstrates proper QC and patient testing procedure of AdenoPlus test <input type="checkbox"/> Correctly interprets patient and quality control results <input type="checkbox"/> Verifies test kit is functioning properly, supplies and reagents are not outdated, and quality controls are within acceptable limits <input type="checkbox"/> Describe Quality Control Remedial Action for out-of-control results			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

Bacterial Vaginosis (OSOM BV Blue) Rapid Test Checklist	# DAHS-NSCBVRT24-POCT		
REFERENCES:		Date	Verifier Initials
1. UC Davis Health Policy POCT.58, Bacterial Vaginosis Rapid Test by OSOM BV Blue Test, Sekisui Diagnostics			
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes proper sample collection of vaginal fluid <input type="checkbox"/> Demonstrates proper test procedure <input type="checkbox"/> Correctly interprets results <input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls). Articulates appropriate action for out-of-control results <input type="checkbox"/> Demonstrates quality control testing using Sekisui OSOM BVBLUE positive and negative sialidase controls. Describes stability of controls and when testing must be performed 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

CoaguChek Checklist	# DAHS-NSCCOAGU24-POCT		
REFERENCES:		Date	Verifier Initials
1. UC Davis Health Policy POCT.42: Prothrombin Time (INR), Whole Blood by Roche CoaguChek XS			
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describe the Roche CoaguChek XS Meter components and their function <input type="checkbox"/> Describe the indications for use <input type="checkbox"/> Describe specimen collection criteria and limitations <input type="checkbox"/> Describe storage requirements, stability criteria, and handling of testing strips <input type="checkbox"/> Demonstrate how to ensure matching of meter and test strips with the code chips <input type="checkbox"/> Describe the integrated quality controls <input type="checkbox"/> Describe remedial action for Out-of-Control Results 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

Cobas® Liat® System Checklist	# DAHS-NSCCOBAS24-POCT		
REFERENCES: 1. UC Davis Health Dept. Pathology & Laboratory Medicine (QC book located in the ED) 2. UC Davis Health Policy POCT.59: Roche Cobas Liat System, Molecular testing of Infectious Disease		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes sample collection and handling <input type="checkbox"/> Describe the stability and handling of the reagent kits <input type="checkbox"/> Verbalizes understanding of valid results. Articulates what to do when an indeterminate or invalid result occurs <input type="checkbox"/> Direct observations of routine patient testing, including, as applicable, patient identification and sample processing and testing <input type="checkbox"/> Demonstrates routine cleaning of instrumentation and immediate work area 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

ColoScreen Checklist	# DAHS-NSCCOLO24-POCT		
REFERENCES: 1. UC Davis Health Policy POCT.19: Occult Blood, Fecal by Helena Laboratories ColoScreen 2. ColoScreen: A Test for Occult Blood, Package Insert, Helena Laboratories		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and associated policy <input type="checkbox"/> Identify factors that may affect the guaiac test <input type="checkbox"/> Describes collection of the stool specimen (able to teach patient for home use) <input type="checkbox"/> Demonstrates proper testing procedure using ColoScreen developer <input type="checkbox"/> Correctly interprets specimen results <input type="checkbox"/> Demonstrates how to perform QC using the internal ColoCheck Monitors. Verbalizes understanding of procedure to follow when control monitors fail <input type="checkbox"/> Describes the stability and handling of the guaiac slides and ColoScreen developer 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:

These skills will be considered complete when all below performance criteria are completed.

Gastrocult Checklist	# DAHS-NSCGAST24-POCT		
REFERENCES:		Date	Verifier Initials
1. UC Davis Health Policy POCT.20: Occult Blood, Gastric Fluid 2. Gastrocult, Package Insert, Beckman Coulter			
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Identify factors that may affect the Gastrocult test <input type="checkbox"/> Demonstrates proper testing procedure of gastric aspirate, using Gastrocult developer <input type="checkbox"/> Correctly interprets results <input type="checkbox"/> Demonstrates how to perform QC using internal Performance Monitors. Verbalizes understanding of procedure to follow when control monitors fail <input type="checkbox"/> Describes the stability and handling of the Gastrocult slides and developer			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:

These skills will be considered complete when all below performance criteria are completed.

HemoCue Hb 201 DM Checklist	# DAHS-NSCHEMOC24-POCT		
REFERENCES: <ol style="list-style-type: none"> UC Davis Health Policy POCT.43: Hemoglobin, Whole Blood by HemoCue Hb 210 DM HemoCue Hb 201 DM Analyzer Instructions for Use Manual, HemoCue HemoCue Hb 201 DM Analyzer & HemoCue DM Docking Station Reference Manual, HemoCue HemoCue Hb 201 Microcuvettes Package Insert, HemoCue HemoCue Cleaner Package Insert, HemoCue HemoCue Technical Letter No 2 September 2010, HemoCueAB R & D Systems R & D Glu/Hgb Control Package Insert, 6/07, R & D Systems, Inc. 		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describe how to handle the microcuvettes and proper sampling technique <input type="checkbox"/> Describe the storage and stability of the liquid controls <input type="checkbox"/> Demonstrate how to perform daily QC using liquid controls <input type="checkbox"/> Describe Quality Control Remedial Action for out-of-control results <input type="checkbox"/> Describe how to perform patient testing on the HemoCue Hb 201 DM Analyzer <input type="checkbox"/> Describe what to do when you have a critical value <input type="checkbox"/> Describe the limitations of the HemoCue HB 201 DM procedure <input type="checkbox"/> Describe the maintenance for the HemoCue Hb 201 DM Analyzer 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

Hemoglobin A1c Afinion™ 2 Analyzer Checklist	# DAHS-NSCHEMOG24-POCT	Date	Verifier Initials
REFERENCES: 1. UC Davis Health POCT Policy 65:Hemoglobin A1c: Whole Blood: Afinion 2 Analyzer			
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy. <input type="checkbox"/> Demonstrates routine patient test performance, including as applicable, patient identification and preparation, specimen collection, handling, processing, and testing. <input type="checkbox"/> Records and reports test results, including reporting critical results as applicable. <input type="checkbox"/> Describes limitations of the procedure <input type="checkbox"/> Performs instrument maintenance and function tests, as applicable. <input type="checkbox"/> Describes problem solving process. 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

Lactate Plus Checklist	# DAHS-NSCLPMSM24-POCT	Date	Verifier Initials
REFERENCES: 1. Department of Pathology & Laboratory Medicine POCT Manual POCT.66			
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and associated policy. <input type="checkbox"/> Direct observations of routine patient test performance, including, as applicable, patient identification and preparation, specimen collection, handling, processing and testing. <input type="checkbox"/> Monitoring the recording and reporting of test results, including as applicable reporting of critical results. <input type="checkbox"/> Describes the limitations of the procedure. <input type="checkbox"/> Direct observation of performance of instrument maintenance and function checks, as applicable. <input type="checkbox"/> Evaluation of problem-solving skills. 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed	

Medtronic HMS Plus 2024	# DAHS-NSCMEDHMSP24-POCT	Date	Verifier Initials
References: 1. UC Davis Health Policy POCT.63:Medtronic HMS Plus			
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy. <input type="checkbox"/> Direct observations of routine patient test performance, including, as applicable, patient identification and preparation, specimen collection, handling, processing, and testing. <input type="checkbox"/> Monitors the recording and reporting of test results, including as applicable reporting of critical results. <input type="checkbox"/> Direct observations of performance of instrument maintenance and function checks as applicable. <input type="checkbox"/> Assessment of test performance through previously analyzed specimens, internal blind testing samples or external proficiency testing samples. <input type="checkbox"/> Evaluation of problem solving skills 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING	
Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed	

Nova StatStrip Checklist	# DAHS-NSCNOVA24-POCT		
References: 2. Instructions for Use: Nova Biomedical StatStrip Glucose Hospital Meter System. Version 1.86. UC Davis Health 3. UC Davis Health Policy POCT.51: Glucose, Whole Blood by Nova StatStrip Glucose Hospital Meter System		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes the Nova StatStrip Glucose Hospital Meter System components and their functions <input type="checkbox"/> Describes the stability and handling of glucose test strips and controls, including open dating, expiration dates and procedure notes <input type="checkbox"/> Describes specimen requirements: amount, types, and recommended anticoagulants <input type="checkbox"/> Describes the conditions when it is not appropriate to perform a glucose fingerstick <input type="checkbox"/> Demonstrates how to perform quality control testing using one control solution. Describes when QC must be run, what QC Lockout is, and what to do if the QC test results FAIL <input type="checkbox"/> Demonstrate how to run a patient test. (Using one control run as a patient) <input type="checkbox"/> Describes the procedure for critical values and questionable patient test results <input type="checkbox"/> Describes the limitations of the Nova StatStrip Glucose Hospital Meter System glucose test results <input type="checkbox"/> Demonstrates how to review results on the meter <input type="checkbox"/> Describes the routine maintenance for the Nova StatStrip Glucose Hospital Meter System, docking station, and accessory box 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING			
Name:		Employee ID#:	
Unit:		Title:	
These skills will be considered complete when all below performance criteria are completed.			
pH Paper Checklist		# DAHS-NSCPHP24-POCT	
REFERENCES: 1. UC Davis Health Policy POCT.22: pH Fluid Nitrazine Paper 2. UC Davis Health Policy POCY.23: pH Fluid by pHdracid Paper 3. Instructions for Use: pHizatest Product Insert. Rev. Nov 15, 2019 4. pHizatest Specification Sheet, version #4, 06-19-2013. 5. Micro Essentials pHizatest Product Information Sheet, Rev. 12-20-12. 6. Micro Essential Hydrion Product Information Sheet, Rev. 10-16-13.		Date	Verifier Initials
<input type="checkbox"/> Completes online module and reads associated policies <input type="checkbox"/> Describes indications of use for each paper type <input type="checkbox"/> Describes the stability and handling, including open dating, expiration dates and quality control <input type="checkbox"/> Describes/demonstrates how to perform and result a test <input type="checkbox"/> Describes the limitations of the procedure			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:
These skills will be considered complete when all below performance criteria are completed.	

Urinalysis dipstick Clinitek Status+ Connect	# DAHS-NSCURID24-POCT		
REFERENCES: <ol style="list-style-type: none"> 1. UC Davis Health Policy: POCT.54: Urine Dipstick by Siemens Uristix and Multistix 10SG Reagent Strips 2. MAS UA Control, Package insert, Thermo Scientific. 3. Multistix 10 SG Reagent Strips, Package insert, Siemens. 4. Uristix Reagent Strips, package insert, Siemens. 5. Keto-Diastix Reagent Strips, Siemens. 		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Collects appropriate urine specimen in clean, dry container <input type="checkbox"/> Demonstrates proper test procedure <input type="checkbox"/> Properly compares the reagent areas to the corresponding color chart on the bottle label at the times specified <input type="checkbox"/> Demonstrates quality control testing using the MAS UA controls. Describes stability of controls and when controls must be run <input type="checkbox"/> Verbalizes understanding of procedure to follow when control results are out-of-range 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING			
Name:		Employee ID#:	
Unit:		Title:	
These skills will be considered complete when all below performance criteria are completed			
Urine Drug Screen by Abbott iCup®		# DAHS-NSCALERE24-POCT	
REFERENCES:		Date	Verifier Initials
1. UC Davis Health Policy POCT.60: Urine Drug Screen by Abbott iCup DX Pro 2			
<input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Describes proper sample collection of urine and demonstrates importance of appropriate collection times <input type="checkbox"/> Demonstrates or describes proper test procedure <input type="checkbox"/> Correctly interprets results <input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls) and articulates appropriate action for out-of-control results <input type="checkbox"/> Describes external quality control testing procedure, stability of reagents and requirements of performance			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years

2024 POINT OF CARE TESTING

Name:	Employee ID#:
Unit:	Title:

These skills will be considered complete when all below performance criteria are completed.

Urine Pregnancy Test Checklist 2023	# DAHS-NSCURIP23-POCT		
REFERENCES: 1. UC Davis Health Policy POCT.41: hCG Pregnancy Test, Urine by Sure-Vue hCG 2. Sure-Vue® Urine hCG, Package Insert, Fisher Scientific Company 3. Sure-Vue® hCG Urine Control Set, Package Insert, Fisher Scientific Company		Date	Verifier Initials
<ul style="list-style-type: none"> <input type="checkbox"/> Completes online module and reads associated policy <input type="checkbox"/> Collects appropriate urine specimen in clean, dry container. Equilibrate urine and test device to room temperature <input type="checkbox"/> Demonstrates proper test procedure <input type="checkbox"/> Correctly interprets results <input type="checkbox"/> Verbalizes understanding of valid results (internal procedural controls). Articulates what to do when there is no control line visible <input type="checkbox"/> Demonstrates quality control testing using Sure-Vue positive and negative hCG controls. Describes stability of controls and when they must be run 			

DO NOT SEND CHECKLIST TO CPPN
File Checklists in Home Department for 2 Years