

**Zoll R Series ALS Skills Checklist #DAHS-NSCRALS17**

<b>Name:</b>	<b>Employee ID #:</b>
<b>Unit:</b>	<b>Title:</b>

PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.  
 These skills will be considered complete when all below performance criteria are completed and pages 1, 2, and 3 have been scanned and emailed to: [cppn@health.ucdavis.edu](mailto:cppn@health.ucdavis.edu)

<b>Method of Instruction</b>	<b>Preceptor Verification of Skill: Method of Evaluation</b>
CP=Clinical Practice D=Demonstration E=Education Session OM=Online Module P=Policy/ Procedure Review SP=Study Packet REF=Reference Information only – not assigned in UCL	N/A=Not applicable to specific patient care area O=Observation (in clinical setting) OT=Online Test RD=Return Demonstration T=Written Test V=Verbal

<b>Zoll R Series ALS Skills Checklist #DAHS-NSCRALS17</b>	<b>Method of Instruction</b>	<b>Method of Evaluation</b>	<b>Date</b>	<b>Verifier Initials</b>
<b>References:</b> 1. <a href="#">UC Davis Health Policy 6005: Automated External Defibrillator (AED-Zoll)</a>				
Completed the assigned ZOLL R Series ALS Defibrillator tutorials in UC Learning				
<b>TEST MODE</b>				
Successfully demonstrates 30 Joule defibrillator test				
Can check and change paper				
<b>AED/MANUAL MODE</b>				
Can turn on device and convert from AED to manual mode				
<b>MONITOR MODE</b>				
Applies 3-lead or 12-lead ECG				
Locates Recorder key and prints a strip				
Access HR menu and demonstrate how to change settings				
Locate NIBP soft key and activate manual BP measurement				
Access NIBP menu and verbalize options				
Demonstrate how to change NIBP alarm settings				
Change NIBP mode from Manual to Automatic				
Change Automatic mode intervals				
Access SpO2 menu and verbalize options				
Access CO2 menu and verbalize options				
<b>MANUAL DEFIBRILLATION</b>				
Locates multifunction cable				
Confirms shockable rhythm				
Selects defibrillator mode (red)				
Presses Energy Select or Charge button				

**Zoll R Series ALS Skills Checklist #DAHS-NSCRALS17**

<b>Name:</b>	<b>Employee ID #:</b>
<b>Unit:</b>	<b>Title:</b>

PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.  
 These skills will be considered complete when all below performance criteria are completed and pages 1, 2, and 3 have been scanned and emailed to: [cppn@health.ucdavis.edu](mailto:cppn@health.ucdavis.edu)

<b>MANUAL DEFIBRILLATION continued</b>	<b>Method of Instruction</b>	<b>Method of Evaluation</b>	<b>Date</b>	<b>Verifier Initials</b>
Tells everyone to stand clear				
Delivers shock at desired energy level				
Defines and adjusts energy levels for adults (120,150, 200J) and Pediatrics (2-4J/kg)				
<b>CPR FEEDBACK</b>				
Demonstrates steps to fill CPR Index™ – understands proper rate/depth				
Shows that if rate is too slow, metronome beeps and Rate prompt appears				
Speeds up to silence metronome and allow the Rate prompt to disappear				
Shows that if depth is too shallow, the Depth prompt appears on the screen				
Pushes hard to allow Depth prompt to disappear				
Demonstrates understanding of See-Thru CPR® filtered ECG				
<b>SYNCHRONIZED CARディオVERSION</b>				
Puts device into SYNC mode				
Selects desired energy				
Presses charge button				
Tells everyone to stand clear				
Delivers synchronized shock				
States and demonstrates that SYNC must be activated for each and every synchronous cardioversion				
<b>PACING</b>				
Turns up pacing output (mA) until capture is achieved – identifies capture				
Adjusts pace rate				
Understands pausing for visualizing patients underlying rhythm				
<b>PADS</b>				
Connects OneStep™ pads to OneStep cable (or other pads/paddles if applicable)				
Opens OneStep packaging correctly				
Demonstrates proper pad placement for defibrillation, pacing, and cardioversion				
Identifies CPR Sensor and explains its purpose				

**Zoll R Series ALS Skills Checklist #DAHS-NSCRALS17**

<b>Name:</b>	<b>Employee ID #:</b>
<b>Unit:</b>	<b>Title:</b>

PERFORMANCE CRITERIA - Unless otherwise specified all skills will be demonstrated in accordance with the appropriate UC Davis Health Policy and Procedure.  
 These skills will be considered complete when all below performance criteria are completed and pages 1, 2, and 3 have been scanned and emailed to: [cppn@health.ucdavis.edu](mailto:cppn@health.ucdavis.edu)

INTERNAL PADDLES	Method of Instruction	Method of Evaluation	Date	Verifier Initials
Understands how to connect internal paddles to OneStep™ cable				
Selects defibrillator mode (red)				
Understands 10J default energy level with range of 1 to 50 Joules				
<b>SUPERUSER/TRAINER</b>				
Demonstrate how to use additional options (Mentor mode, Set the clock, etc.)				
Understands how to change parameter settings (NIBP, EtCO2, SpO2)				
Understands purpose of Code Marker				
Can access data from the code (Print Chart, Print Log, or Transfer Data).				
User demonstrates sufficient understanding of device to train other users in its use.				

**PRECEPTOR SIGNATURE:**

Signature and Printed Name of Verifier (preceptor or other verified personnel) who have initialed on this form:

Initials:	Print Name:	Signature:

**PRECEPTEE STATEMENT AND SIGNATURE:**

I have read and understand the appropriate UC Davis Health Policies/Procedures and/or equipment operations manual, I have demonstrated the ability to perform the verified skills as noted, and I have the knowledge of the resources available to answer questions.

<b>Name:</b>	<b>Signature:</b>	<b>Date:</b>
--------------	-------------------	--------------